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
# Workflow broken Be-window

**Call the DESY Technical Emergency Service (2500, 5555) and the FS Shift Service (3868)**

Call Technical hall service and FS-Safety, outside of office hours per e-mail to [photon-science.safety-managers@desy.de](mailto:photon-science.safety-managers@desy.de) and Fselog

- **First safety measure: block area, enter with PPE set only:** overshoes, mouth protection, overall, gloves. Use of an adhesive mat at the lab entrance
- **On site: risk analysis by TN + User/Staff..**
  - Be-window inside detector
  - Be-window under vacuum
  - Open Be-window at experiment
- **Safety measures depending on risk and competence of staff on site, e.g.:**
  - Keep the area blocked, if necessary request wipe test or decontamination via V8 safety department.
  - Pick up visible splinters with adhesive tape
  - Check splinters for completeness
  - Use of special beryllium vacuum cleaner (ask shift service/ technical hall service for Be vacuum cleaner, adhere to operating instruction!)
  - Tape up the windows, pack airtight.
  - Dispose of splinters, protective clothing and adhesive mat as waste, if necessary.
  - Send in detector to company
  - Check whether future substitution of beryllium is possible or a technical solution to prevent a crash.
  - ...

Information to	Phone, E-mail
Technical emergency service (TN) Shift service	2500, 5555 3868
FS-Safety Manager	<a href="mailto:photon-science.safety-managers@desy.de">photon-science.safety-managers@desy.de</a>
FS-eelogbook	<a href="http://tfinfo.desy.de/FSelog/">http://tfinfo.desy.de/FSelog/</a> Done by shift service
Wipe test, decontaminaton	via V8 safety department